

CONSUMABLE & SYSTEM ITEMS

Catalogue



Consumables

XGD TIP HOLDER

XGD General Tip holder



For use on XGD probes with R-Type Tips

GTH-25

QTY :25

R-TYPE PROBE EAR TIPS (for XGD and UGD probes)

DP + TE

R4.8S



Fits - 4.8mm ear canal
Small newborns

T-R4.8S

QTY: 100

R5.8B



Fits - 5.8mm ear canal
Most newborns

T-R5.8B

QTY: 100

R6.8B



Fits - 6.8mm ear canal
Large newborns and first year infants

T-R6.8B

QTY: 100

R7M



Fits - 7mm ear canal
Infants and children

T-R7M

QTY: 100

R8M



Fits - 8mm ear canal
Infants and small adult ears

T-R8M

QTY: 100

DPC COUPLER TUBES

UGD Probe Couplers



For use on UGD probes with R-Type Tips

QTY : 10

QTY : 100

DPC-10

DPC-100

R9M



Fits - 9mm ear canal
Most adult ears

T-R9M

QTY: 100

R11M



Fits - 11mm ear canal
Large adult ears

T-R11M

QTY: 100

R13M



Fits - 13mm ear canal
Extra large adult ears

T-R13M

QTY: 100

Small Foam Tip



Small - for pediatric and adult ears

T-R100S

QTY: 100

Large Foam Tip



Large - for adult ears

T-R100L

QTY: 100

TPC COUPLER TUBES

UGS Probe Couplers



For use on UGS probes with T-Type Tips

QTY : 10

QTY : 100

TPC-10

TPC-100

T3E



Fits - 3mm ear canal
Small and premature newborns

T-T3E

QTY: 100

T4.5C



Fits - 4.5mm ear canal
Small newborns

T-T4.5C

QTY: 100

T5C



Fits - 5mm ear canal
Newborns

T-T5C

QTY: 100

T5.5B



Fits - 5.5mm ear canal
Most newborns

T-T5.5B

QTY: 100

T6.5B



Fits - 6.5mm ear canal
Large newborns and first year infants

T-T6.5B

QTY: 100

T7M



Fits - 7mm ear canal
Infants and children

T-T7M

QTY: 100

T8M



Fits - 8mm ear canal
Infants and small adult ears

T-T8M

QTY: 100

T9M



Fits - 9mm ear canal
Most adult ears

T-T9M

QTY: 100

T11M



Fits - 11mm ear canal
Large adult ears

T-T11M

QTY: 100

T13M



Fits - 13mm ear canal
Extra large adult ears

T-T13M

QTY: 100

Large Foam Tip



Large - for adult ears

T-T200L

QTY: 100

ABR CONSUMABLES

Electrodes



Snap
Tab

ABR-E-SNAP
ABR-E-TAB

QTY:25
QTY:60

Ear Cups



ABR-CUP

QTY: 20

MAXIMISE PROBE LIFE

It is good practice to change the **TIP HOLDER OR COUPLER** regularly, for example, after every **10 EAR TIP** changes. The **TIP HOLDER OR COUPLER** should be replaced immediately if seen to be contaminated. Always ensure that your local infection control policy is followed.

Otodynamics probe **EAR TIPS** are disposable and **MUST** be discarded after each patient.



System Items

XGD PROBE

XGD DP+TEOAE



For use with Otoport / Otocheck
ABR TE+DP testing

1 metre PR-XGD QTY : 1

XGD Left/Right Identification Labels



For use with Otoport / Otocheck.
ABR binaural left & right identification
for ABR + DP+TE testing

XGD-RL QTY : 2

PROBES

UGD DP+TEOAE



For use with Echoports
and Otoport / Otochecks

2 metre PR-UGD QTY : 1

UGD DP+TEOAE (red)



For use with Otoport / Otocheck ABR binaural
right identification for ABR+DP+TE testing

1 metre PR-ABRRD QTY : 1

UGD DP+TEOAE (blue)



For use with Otoport / Otocheck ABR binaural
left identification for ABR+DP+TE testing

1 metre PR-ABRLD QTY : 1

UGS TEOAE



For use with Echoports
and Otoport / Otochecks

2 metre PR-UGS QTY : 1

UGS TEOAE (red)



For use with Otoport / Otocheck ABR binaural
right identification for ABR +TE testing

1 metre PR-ABRRS QTY : 1

UGS TEOAE (blue)



For use with Otoport / Otocheck ABR binaural
left identification for ABR+TE testing

1 metre PR-ABRLS QTY : 1

PROBE BODY & LID

TE Claret Body & Lid



For UGS probes

BGS QTY : 10

TE Red Body & Lid



For UGS probes
(Right identification)

BGS-R QTY : 10

TE Blue Body & Lid



For UGS probes
(Left identification)

BGS-L QTY : 10

DP Grey Body & Lid



For UGD probes

BGD QTY : 10

DP Red Body & Lid



For UGD probes
(Right identification)

BGD-R QTY : 10

DP Blue Body & Lid



For UGD probes
(Left identification)

BGD-L QTY : 10

PROBE ITEMS

Probe Pouch



PR-POUCH QTY : 10

Probe Cable Clip



PR-CLIP QTY : 10

1cc OAE+ABR Test Cavity



ABR-CAV QTY : 1

1cc OAE+ABR Test Cavity use with crocodile clip electrodes



NHSP-ABR-CAV QTY : 1

OTOPORT ABR STAND

ABR Desktop Stand/Crib Hook



Desktop stand or hooking
to crib for Otoport ABR

ABR-DS QTY : 1

OTOPORT OAE STAND

Docking Station



(Please Note: This Docking Station and Charging Cradle
are not compatible with the OAE+ABR instrument)

ODS QTY : 1

Charging Cradle



OCC QTY : 1

INFECTION CONTROL

Otoport ABR Infection Control Sleeve



ABR-INF QTY : 10

Otoport OAE Infection Control Sleeve



OP-INF QTY : 10

System Items

USB CABLES

USB Cable



0.4 metre USB-CAB-04 QTY : 1
2 metre USB-CAB-2 QTY : 1

OTOPORT SYSTEM ITEMS

Download Cable



Connecting
Otoport to PC

OP-CAB QTY : 1

Charger



Compact wired
Otoport charger
(including mains cable)

OP-CHG QTY : 1

ABR CABLES

Stud Electrode Cable



1 metre ABR-EC-1 QTY : 1
2 metre ABR-EC-2 QTY : 1

Crocodile Clip Electrode Cable



1 metre ABR-EC1-T QTY : 1
2 metre ABR-EC2-T QTY : 1

OTOPORT MINI PRINTER & ITEMS

Mini-Printer



Including Charger, Printer
Cable and Paper

PP QTY : 1

Mini-Printer Paper Rolls



Standard printer
paper rolls

OMP-PAP QTY : 10

Mini-Printer Cable



Connecting Otoport
to Printer

OMP-CAB QTY : 1

Mini-Printer Charger



Replacement
printer charger

OMP-CHG QTY : 1

Mini-Printer Battery Pack



Replacement
battery pack

OMP-BATT QTY : 1

OTOPORT FABRIC CARRY CASE

Otoport OAE (standard)



OP-CAS QTY : 1

Otoport OAE with printer compartment



OPP-CAS QTY : 1

Otoport ABR



ABR-CAS-ST QTY : 1

Otoport ABR with printer compartment



ABR-CAS-LG QTY : 1

Software

ILOV6



ILO v6 for Otoport and Echoport
Advanced Data Management and
Analysis Software - Single
Instrument Licence
(includes EZ•Screen)
Compatible with Windows 11 and 10

ILO V6

EZ•SCREEN



EZ•Screen for Echoport
Data Management
Single Instrument Licence
Compatible with
Windows 11 and 10

EZ-SCREEN

OTOLINK



Manage your data easily and
efficiently with the Otolink
PC patient data management suite
Single Instrument Licence
Compatible with Windows 11 and 10

OTOLINK

NOAH



ADD NOAH
Our NOAH option makes Otolink
PC data management suite fully
compatible with Noah 4

NOAH



Otodynamics Ltd.

30-38 Beaconsfield Road Hatfield Herts AL10 8BB UK

T. +44 (0)1707 267540 E. info@otodynamics.com

www.otodynamics.com